

Alexa Fluor 647 Mouse IgG2b 同型对照

目录号	规格
FMCK003-01-025	25 次
FMCK003-01-100	100 次

克隆号: 4A2B5.3

亚型: Mouse IgG2b

缓冲液: 含有 0.2%BSA 和 0.1%Proclin 300 磷酸盐缓冲溶液, PH7.2-7.4

应用: 流式细胞术

推荐用量: 5 μ l/次

激光器: 红色 (633 nm)

最大吸收/发射波长: 633nm/668nm

存储条件: 4 度, 避光保存

简介:

本抗体亚型是 mouse IgG2b。可作为流式实验中同型对照使用。

实验方法:

流式细胞术表面染色方案

流式细胞术胞内染色方案

相关产品:

Cat NO	Name
FXP001	RBC Lysing Buffer (10x)
FXP005	Cell Staining Buffer
FXP006	BFA, Brefeldin A
FXP008	Fixation Buffer
FXP009	Permeabilization Wash Buffer
FXP045	Monensin Solution (1,000X)
FXP048	Intracellular Cell Stimulation Cocktail (500X) 4 x 100 uL
FXP049	Intracellular Cell Stimulation Cocktail (plus protein transport inhibitors) (500X) 100 uL
FXP055	Sample preservation buffer
FXP064	Human Fc Receptor Blocking Solution
FMU16/32	Purified anti-mouse CD16/32

相关文献:

1. Niemann J, et al. 2019. Nat Commun. 10:3236.
2. Schroder WA, et al. 2019. Sci Rep. 9:12421.
3. Deniset JF et al. 2019. Immunity. 51(1):131-140 .
4. Hong Y, et al. 2019. J Extracell Vesicles. 8:1670893.
5. Yang C, et al. 2019. J Transl Med. 1.004166667.
6. Miao L, et al. 2020. Theranostics. 0.7625.
7. Zhang M, et al. 2019. Autophagy. 1.423611111.

4A Biotech Co., Ltd.

Add: No.Room 201, building 17, No. 8, Jinfeng Road, science and Technology City, Gaoxin District, Suzhou
www.4abio.com , Technical support: tech@4abio.com