

APC 标记山羊抗人 IgG (MSA)

目录号 规格

SPA1031-100 100ug

SPA1031-1000 1mg

抗体类型 多克隆

缓冲液: 含有 0.2%BSA 和 0.1%Proclin 300 磷酸盐缓冲溶液, PH7.2-7.4

应用: 流式细胞术

浓度: 1mg/ml

激光器: 红色激光 (633 nm)

最大吸收/发射波长: 650nm/660nm

存储条件: 4 度, 避光保存

简介:

本品是 APC 标记山羊抗人 IgG, 能与兔 IgG(H+L)发生特异性反应, 用于流式实验。

实验方法:

1.取 EDTA 抗凝的外周血 100 μ l, 加入到 5ml 流式管管底;

4A Biotech Co., Ltd.

Add: No.Room 201, building 17, No. 8, Jinfeng Road, science and Technology City, Gaoxin District, Suzhou
Website: www.4abio.net Toll Free (China only): 400-7060-959 Technical support: tech@4abio.com

- 2.加入适量纯化抗体，与外周血混匀后，室温孵育 20 分钟；
- 3.加入 2ml 1×红细胞裂解缓冲液，与样本混合后孵育 10 分钟，裂解红细胞（推荐 cat: FXP001）；
- 4.300g 离心 5 分钟，弃上清；
- 5.加入 2ml 1×PBS 缓冲液重悬细胞，然后 300g 离心 5 分钟，弃上清，用 100μl PBS 重悬细胞；
- 6.加入适量 APC 标记山羊抗人 IgG，室温避光孵育 20 分钟；
- 7.重复步骤 5 一次；
- 8.加入 0.5ml 1×PBS 缓冲液重悬细胞，通过流式细胞仪进行检测进行
- 注：**未避免荧光猝灭，染色完成的样本应在 2 小时内上机测定，若不能及时上机检测，可将样本经细胞固定液固定后，置于 2~8℃条件下避光保存，并于 24 小时内完成上机检测。

相关产品:

目录号	名称
25610	Lymphocyte separation solution
FXP001	RBC Lysing Buffer (10x)
FXP005	Cell Staining Buffer
FXP055	Sample preservation buffer
FXP064	Human Fc Receptor Blocking Solution

4A Biotech Co., Ltd.

Add: No.Room 201, building 17, No. 8, Jinfeng Road, science and Technology City, Gaoxin District, Suzhou
Website: www.4abio.net Toll Free (China only): 400-7060-959 Technical support: tech@4abio.com